

Mistral Multi-Mixers®

- ♦ Eleven specialized, interchangeable attachments allow for the performance of a wide variety of laboratory mixing tasks
- ♦ Compact size conserves valuable bench top space
- ♦ Speed is adjustable from 50 to 2400 rpm with 0.120" (0.3 cm) circular motion
- ♦ Heavy duty DC motor
- ♦ Counterbalanced drive mechanism with flexible suspension—mounted on rubber feet to increase stability
- ♦ Rugged steel body

construction with baked enamel finish

- ♦ Mixer attachments are constructed of corrosionresistant anodized aluminum and can be changed in seconds
- ♦ CE and CSA certified

| ORDERING INFORMATION | | | | | |
|------------------------------|--------|-------|--|----------------------|-----------------------|
| Model | Volts* | Watts | Dimensions W" x L" x H" (cm) | Net Wt. Ibs. (kg) | Ship Wt. Ibs. (kg) |
| 4600 | 120 | 120 | 81/8 x 101/2 x 7 (20.6 x 26.7 x 17.8) | 18 (8.2) | 22 (10) |
| 4600-1 | 240 | | | | |
| *Both models rated 50/60 Hz. | | | | | |

| ACCESSORIE | 3 | | | |
|------------|---|--|--|--|
| Model | Description | | | |
| 4601 | Utility Tray, $6\frac{1}{2}$ "L x 12"W (16.5 x 30.5 cm) with rubber nonslip mat | | | |
| 4602 | Universal Holder Attachment for flask, beakers and bottles up to 500 ml capacity | | | |
| 4603 | Separatory Funnel Attachment for 1 funnel up to 500 ml capacity | | | |
| 4604 | Adjustable Bottle Holder Attachment, horizontally holds bottles up to 500 ml capacity | | | |
| 4605 | Four-Place Holder Attachment, tri-grip clamps hold round bottom flasks or bottles up to 250 ml capacity | | | |
| 4606 | Test Tube Holder Attachment holds 36 13 mm test tubes | | | |
| 4607 | Round Bottom Flask Holder Attachment, vertically holds 1 flask or bottle up to 500 ml capacity | | | |
| 4608 | Rubber Pad Attachment for hand-held test tubes 45 %" x $5\frac{1}{2}$ " (11.7 x 14 cm) | | | |
| 4609 | Holder Attachment holds Model 4610 Sieve Set | | | |
| 4609-5 | Holder Attachment holds Model 4610-5 Sieve Set | | | |
| 4610 | Sieve Set, 2" (5 cm) diameter, 4 stackable holders, #25, 35, 45, 60, 80, 120, 170 and 230 mesh screens | | | |
| 4610-5 | Sieve Set, 3" (7.6 cm) diameter, 4 stackable holders, #18, 35, 60 and 120 mesh screens | | | |
| 4611 | Microtube platform holds 64 0.5 ml and 81 1.5 ml microtubes | | | |
| 4612 | Replacement Foam Pad for 4611 platform | | | |